

Patient Case History

- 48 year old, Female Patient, attending hospice for support with breast cancer, for which she has had surgery.
- Patient has Lymphoedema and needs support with lymphatic drainage.
- The patient is self-caring and has previously not worn a compression bra for her condition

Product and use

The Patient wore the Solidea Silva wave compression bra*. Which is a medical-grade compression garment designed to alleviate and manage lymphatic circulation disorders, particularly lymphoedema in the breast and armpit area. It has a micro-massaging fabric reactivates circulation and helps eliminate fluids

The patient has been wearing the bra daily (not overnight) for 5 months. The Patient wore a XXL bra.

Patient and Clinician were asked to provide feedback on the bra after use.



Patient Feedback

(1 = Poor - 5 = Excellent)

How would you rate the overall fit of the garment?	4
What would the patient rate the breathability of the fabric?	5
How well did the product manage moisture during wear? (e.g. sweat absorption):	5
How would you rate the durability of the product? (e.g. quality of materials, resistance to tearing)	5
Did the garment stay in place while in use?	5
What was the patients' level of comfort while wearing the garment?	5
How would they rate the texture and feel of the garment?	5
How easy did the patient find of application?	5
How did the patient find the ease of removal?	5
How would you rate your mobility while using the garment?	5

Clinician Feedback

(1 = I strongly disagree - 10 = I agree)

The compression level of the garment was sufficient for the condition.	8
The garment was easy to apply	8
The garment was easy to remove	8
The garment was beneficial/effective for my patient	8
The patient was able to self-care with this garment	8
I would use the product again	8
I was satisfied with the product	8
There was an improvement in the patients condition	8
There was an improvement in the patient's quality of life	9
There was an improvement in the patients skin (if applicable)	9

Pain:

- The patient did not report any changes in pain.

Other Outcomes:

- The quality of the patient's skin improved over the duration of wear.
- Mobility/range of motion improved.
- Bra gave good support

SOLIDEA®
**SILVER
WAVE**

Summary:

The Patient and the clinician were very satisfied with the product and have continued to use this garment to support their care. This care report demonstrates the benefits of this product in managing lymphoedema in the breast. The clinician also noted an improvement to the patient's quality of life and that the patient was able to use this daily and continue to self-care.