Daylong

Need help with measuring?

For more

information please scan here

Credalast Velvet and Classic Made to Measure form

| PART 1 | | |
|--|-------------------------------------|--|
| Name of pharmacy: | | |
| Pharmacy contact name: | | |
| Pharmacy address: | | |
| | | |
| | | |
| Postcode: | | |
| Pharmacy email: | | |
| Pharmacy contact number: | | |
| Patients name: | | |
| Customer account number: | | |
| Previous ref No. (if repeat): | | |
| Circumference Measurements CG CF | Length Measurements | |
| CA measurement must be completed for closed and open toe | | |
| PART 2 Please complete ALL se | ctions and tick the appropriate box | |
| Section One – Choose which range of stocking: | | |
| Credalast Classic | | |
| Credalast Velvet | | |
| | | |
| Section Two – Compression: | Section Four – Foot: | |
| Class I (14-17mmHg) | | |
| Class II (18-24mmHg) | Closed Toe | |
| Class III (25-35mmHg) | | |
| Section Three – Colour: | Section Five – Style: | |
| Natural Below Knee | | |
| Black | Thigh Length | |

| PART 3 Must be completed with numbers or zero in each box | | | | |
|---|--|-----------|------------|--|
| Quantity: | Left | Leg | Right Leg | |
| PART 4 CA measurement must be completed for closed and open toe | | | | |
| Thigh Length (Please complete all below measurements) | | | | |
| C = CIRCU | JMFERENCE | Left (cm) | Right (cm) | |
| CG (Thig | lh) | | | |
| CF (Mid | Thigh) | | | |
| CD (Belo | ow Knee) | | | |
| CC (Calf |) | | | |
| CB (Ank | le) | | | |
| | mference nd the heel) | | | |
| | mference of foot ured at base of toes) | | | |
| L = LENG | ΓH | Left (cm) | Right (cm) | |
| LG (Thig | lh) | | | |
| LF (Mid | Thigh) | | | |
| LD (Belo | ow Knee) | | | |
| LC (Calf |) | | | |
| LB (Ank | le) | | | |
| | en toe) | | | |
| LZ (Clos | sed toe) | | | |
| Below Knee (Please complete all below measurements) | | | | |
| CA measurement must be completed for closed and open toe | | | | |
| C = CIRCU | JMFERENCE | Left (cm) | Right (cm) | |
| CD (Belo | ow Knee) | | | |
| CC (Calf |) | | | |
| CB (Ank | le) | | | |
| | mference nd the heel) | | | |
| | mference of foot sured at base of toes) | | | |
| L = LENG | | Left (cm) | Right (cm) | |
| | ow Knee) | | | |
| LC (Calf | | | | |
| LB (Ank | | | | |
| · · · | en toe) | | | |
| LZ (Clos | sed toe) | | | |
| Additional Notes: | | | | |

Submission of this form confirms that the patient's consent has been obtained to submit these details to the manufacturer, in order to produce the patient's made to measure stockings. * Items will be dispatched within 5 working days if all details are correct and subject to factor's beyond our control Exclusively available from



Guidance notes for Credalast Velvet and Credalast Classic Made to Measure Stockings

Credalast Circular Knit Made to Measure Compression Stockings. Available on Drug Tariff

Credalast stockings are available on **Drug Tariff FP10**

Please ensure **Parts 1-4** are completed overleaf along with the prescription (each option below must be listed on the prescription)

- Velvet or Classic Compression Class Colour Open or Closed Toe
- Below Knee or Thigh Length Quantity for each leg

If you are a **Pharmacist**: email measurement form only

to mtm@daylongdirect.co.uk (scans or smart phone pictures accepted)

If you are a **Patient** or **Nurse Prescriber:** send measurement form and prescription to

FREEPOST RTHL-YUAG-GCST Daylong Direct, Credenhill Ltd, Unit 5 Soloman Road, Cossall Industrial Estate, Ilkeston, Derbyshire, DE7 5UD

or send the prescription electronically to us as the appliance contractor using the NACS code FJ708 and email the form to mtm@daylongdirect.co.uk

Stockings will be dispatched within **5 working days** of Daylong Direct receiving the form, providing there are no queries. If you need any further help please ring us on 0115 932 0144, our experienced team will be able to help.

Daylong Direct currently offers a choice of two stockings, Credalast[®] Velvet and Credalast[®] Classic both developed from premium Swiss stocking manufacturing technology.

Credalast[®] Velvet

Credalast Velvet

class 2, natural, open toe, Below

MTM


NHS

The first prescription Made to Measure stocking available with the Tactel® Climate Effect which wicks away moisture and improves air permeability, helping to keep an optimum skin balance regardless of activity.

Credalast® Classic

This traditional stocking offers a more conventional stiffness and style which suits many patients.

Ŏ The Tactel" The mate effect 00 0 Tacteľ 000



Application Aids

Some users may find it difficult to apply compression hosiery. For this reason we stock a range of highly useful application aids such as the Magnide (open and closed toe stockings) and Easy Slide (open toe only) to help improve the wearer experience.

Credalast® Classic and Credalast® Velvet are registered trade marks. Tactel is a registered trade mark of Invista. © 2024. Credenhill. All Rights Reserved.

Accessories

Thigh length stockings will need suspending as they are not manufactured with a self supporting top band.

We offer the Credalast Suspender Belt which is available on Drug Tariff FP10 or the Credalast Y Shaped Suspenders. We can also provide Coltex - A skin adhesive which is designed to hold compression hosiery in place which can also be ordered via Daylong Direct.

Daylong Direct is a trading name of Credenhill Limited Daylong Direct, Unit 5 Soloman Road,

Cossall Industrial Estate, Ilkeston, Derbyshire, DE7 5UD.