TROM PRIL

# SEALCUFF

IDEAL FOR WET CONDITIONS

**BOTH OF THESE ARE IDEAL FOR WET CONDITIONS...** 



...ONLY ONE IS AVAILABLE ON PRESCRIPTION.

Protect wounds, dressings, PICC lines, bandages, casts and hosiery with **SealCuff® - Complete** waterproof protection, reusable and latex-free.

Exclusively available from





Whether you are managing a patient with a venous leg ulcer, chronic oedema, cancer undergoing intravenous therapy via PICC line, or a child with a broken limb in a plaster cast. SealCuff® is your go-to solution for waterproof protection.

**Available on FP10 prescription**, it ensures safe showering without compromising dressings, bandages, or treatment sites.

#### **Ideal for Use With:**

Wounds & Leg Ulcers - Safeguard advanced dressings during healing

Compression Bandages – Support oedema and lymphoedema patients with bulky wraps

Plaster Casts - Keep casts dry, prevent odour and skin breakdown

**PICC Lines & IV Sites** – Protect insertion areas from water ingress

Compression Hosiery – Maintain garment integrity between changes

#### **Clinician Benefits**

#### **Maintain Clinical Outcomes**

- Reduces risk of avoidable infection by keeping wounds and invasive sites dry
- Minimises dressing changes, conserving time and resources

#### **Enhance Patient Compliance**

- Comfortable neoprene seal promotes regular hygiene routines, reducing non-compliance and associated complications
- Latex-free and skin-friendly, suitable for sensitive skin or long-term use

#### **Deliver Cost-Efficient Care**

- Up to 22-40%\* cost savings vs. the leading alternative
- Durable and **reusable**, lowering long-term costs per patient

# SEALCUFF®

### IDEAL FOR WET CONDITIONS

- Waterproof protection for dressings, PICC lines, casts, hosiery, and wraps.
- Reusable, latex-free, with a gentle yet strong neoprene seal for improved patient comfort and compliance.

### SEALCUFF® ADULT FOOT


### **Limb fit location:** Above ankle

Circumference of limb: (minimum to maximum)

27 to 54cm

#### Length measurements:

Measure from floor to above ankle - max length 30cm





### SEALCUFF® ADULT SHORT LEG

#### Limb fit location:

Just above knee

Circumference of limb: (minimum to maximum) 27 to 54cm

**Length measurements:** 

Measure from floor to just above knee - max length 60cm





### SEALCUFF® ADULT LONG/XL LEG

#### Limb fit location:

Mid thigh

Circumference of limb: (minimum to maximum)

27 to 54cm/32 to 62cm

#### **Length measurements:**

Measure from floor up to mid thigh - max length 70cm Long/XL - max length 100cm





### SEALCUFF® ADULT HAND

Limb fit location:

Mid forearm

Circumference of limb: (minimum to maximum)

16 to 38cm

Length measurements:

Measure from finger tips to mid forearm - max length 35cm





### SEALCUFF® ADULT SHORT ARM

**Limb fit location:** 

Mid bicep

Circumference of limb: (minimum to maximum)

16 to 38cm

**Length measurements:** 

Measure from finger tips to mid bicep - max length 50cm





### SEALCUFF® ADULT LONG ARM

**Limb fit location:** 

Mid to upper bicep

Circumference of limb: (minimum to maximum)

27 to 54cm

**Lenath measurements:** 

Measure from finger tips to mid/ upper bicep - max length 60cm





### SEALCUFF® **PICC LINE**

**Limb fit location:** 

Below and above elbow

Circumference of limb:


(minimum to maximum)

27 to 54cm

**Length measurements:** 

Measure from below and above elbow - max length 40cm





## SEALCUFF®

# Complete Waterproof Protection for Complex Care Needs

#### NHS SUPPLY CHAIN ORDER CODES AND FP10 PRESCRIPTION PIP CODES

SEALCUFF® Waterproof Cast Cover	Circumference of limb (minimum to maximum)	Length measurements	Product Code	PIP Code	NHSSC NPC	
Hand and Arm						
Adult Hand	16 to 38cm	Measure from finger tips to mid forearm - max length 35cm	230131	435-1029	EAZ85093	
Adult Short Arm	16 to 38cm	Measure from finger tips to mid bicep - max length 50cm	230231	435-1094	EAZ85092	
Adult Wide Short Arm	27 to 54cm	Measure from fingertips to mid bicep - max length 60cm	230232	435-1102	EAZ85091	
Adult Long Arm	27 to 54cm	Measure from finger tips to mid/ upper bicep - max length 60cm	230331	435-1110	EAZ85090	
PICC						
IV&PICC Lines-Small	16 to 38cm	Measure from below and above elbow - max length 30cm	231030	435-1177	EAZ85095	
IV&PICC Lines-Large	27 to 54cm	Measure from below and above elbow - max length 40cm	231032	435-1169	EAZ85085	
Foot						
Adult Foot	27 to 54cm	Measure from floor to above ankle - max length 30cm	230631	435-1128	EAZ85083	
Leg						
Adult Short Leg	27 to 54cm	Measure from floor to just above knee - max length 60cm	230531	435-1136	EAZ85088	
Adult Wide Short Leg*	27 to 54cm	Measure from floor to just above knee - max length 60cm	230532	435-1144	EAZ85081	
Adult Long Leg	27 to 54cm	Measure from floor up to mid thigh - max length 70cm	230431	435-1151	EAZ85089	
Adult Long Leg XL	32 to 61cm	Measure from floor up to upper-thigh - max length 100cm	230433	Coming soon	Coming soon	
Child Arm						
Kids Arm - Small	16 to 38cm	Measure from finger tips to mid/ upper bicep - max length 34cm	220230	435-1037	EAZ85087	
Kids Arm - Medium	16 to 38cm	Measure from finger tips to mid/ upper bicep - max length 48cm	220231	435-1045	EAZ85080	
Kids Arm - Large	16 to 38cm	Measure from finger tips to elbow - max length 54cm - elbow to finish - max length 37cm	220232	435-1052	EAZ85086	
Child Leg						
Kids Leg - Small	27 to 54cm	Measure from floor to just above knee - max length 32cm/ Foot max length 27cm	220530	435-1060	EAZ85084	
Kids Leg - Medium	27 to 54cm	Measure from floor to just above knee - max length 48cm / Foot max length 34cm	220531	435-1078	EAZ85082	
Kids Leg - Large	27 to 54cm	Measure from floor up to mid thigh - max length 72cm / Foot max length 34cm	220532	435-1086	EAZ85094	

### SEALCUFF®

## **Complete Waterproof Protection** for Complex Care Needs

Feature	Clinical Benefit		
Complete Waterproof Seal	Secure protection for therapy continuity		
Soft Neoprene Cuff	Flexible and comfortable, ideal for regular use including elderly and frail patients		
✓ Wide Range of Sizes	Suitable for all ages, limb sizes, and complex bandaging needs		
FP10 & NHS Supply Chain Available	Easy to prescribe and reorder		

### **SEAL** IN THE SAVING

28%\* less vs. the main competitor.











Your patient can purchase SealCuff from our website **www.daylong.co.uk** or buy over the phone





Scan the QR Code to order SealCuff®



To find out more visit www.daylongdirect.co.uk

Tel 0115 932 0144 • Email sales@daylongdirect.co.uk