

Solidea® BE YOU Tonic Range – Prescriber Information



The Solidea BEYOU Tonic range combines graduated compression 15–21 mmHg with patented 3D wave-knit micromassage technology. The garments are designed to support venous and lymphatic return, help stimulate microcirculation, and provide a comfortable, wearable option for daily self-management for patients with lipoedema, lymphoedema or chronic oedema.

The range includes 6 options:

BE YOU Tonic, BE YOU Tonic Curvy, BE YOU Tonic Pro, BE YOU Tonic Curvy Pro, BE YOU Capri and BE YOU Capri Curvy.

How to prescribe/order

When prescribing or ordering, specify the full product name, size, colour, PIP code and product code.

Example

Product: Solidea BEYOU Tonic Curvy

Size: 3ML-XL

Colour: Blu Tonic

PIP code: 441-0551

Product code: 0603B5SOSM2603

FEEL THE WAVE

3D Patented Wave-Knit Technology

Micromassage compression




Measuring advice


Measuring and size selection

- 1 Use the relevant size chart for the selected garment type.
- 2 Measure the patient's current height and weight.
- 3 Plot height on the vertical axis and weight on the horizontal axis.
- 4 Select the size where the two measurements intersect.
- 5 Consider the Curvy option where the patient has fuller hips, a curvier body shape, or where standard fit may roll, dig in or be poorly tolerated.
- 6 The garment should feel supportive but not restrictive and should sit flat against the skin.
- 7 Manufacture recommends sizing up if the patient is between sizes, use clinical judgement and prioritise comfort, tolerance, body shape and expected concordance.

BE YOU Tonic Measurement Chart



BE YOU Tonic Curvy Measurement Chart



Exclusively available from



To find out more visit www.daylongdirect.com

For further support please contact us on:

Telephone 0115 932 0144

Email sales@daylongdirect.co.uk

Please scan the QR code to shop online



Solidea® BE YOU Tonic Range Order Codes

Suitable for:

For adult patients with lower-limb lipoedema, lipo-lymphoedema or oedema who need a comfortable, wearable compression option to support daily self-management.

BE YOU Tonic - 15-21 mmHg								
	Blu Tonic		Champagne		Nero		Ribes	
Size	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code
S	0602B5SOSM2601	441-0122	0602B5SOSM2201	441-0171	0602B5SOSM0901	441-0270	0602B5SOSM4601	441-0320
M	0602B5SOSM2602	441-0106	0602B5SOSM2202	441-0155	0602B5SOSM0902	441-0254	0602B5SOSM4602	441-0304
M/L	0602B5SOSM2603	441-0114	0602B5SOSM2203	441-0163	0602B5SOSM0903	441-0262	0602B5SOSM4603	441-0312
L	0602B5SOSM2604	441-0098	0602B5SOSM2204	441-0148	0602B5SOSM0904	441-0247	0602B5SOSM4604	441-0296
XL	0602B5SOSM2605	441-0130	0602B5SOSM2205	441-0189	0602B5SOSM0905	441-0288	0602B5SOSM4605	441-0338

CURVY RANGE – FOR SHAPE DIFFERENTIATION

BE YOU Tonic Curvy - 15-21mmhg								
	Blu Tonic		Champagne		Nero		Ribes	
Size	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code
1 S-XL	0603B5SOSM2601	441-0536	0603B5SOSM2201	441-0577	0603B5SOSM0901	441-0650	0603B5SOSM4601	441-0692
2 M-XL	0603B5SOSM2602	441-0544	0603B5SOSM2202	441-0585	0603B5SOSM0902	441-0668	0603B5SOSM4602	441-0700
3 ML-XL	0603B5SOSM2603	441-0551	0603B5SOSM2203	441-0593	0603B5SOSM0903	441-0676	0603B5SOSM4603	441-0718
4 L-XL	0603B5SOSM2604	441-0569	0603B5SOSM2204	441-0601	0603B5SOSM0904	441-0684	0603B5SOSM4604	441-0726

PRO RANGE – OPEN TOE WITH FOOT

BE YOU Tonic Pro - 15-21mmHg						
	Blu Tonic		Nero		Ribes	
Size	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code
S	0615B5SOSM2601	441-0379	0615B5SOSM0901	441-0429	0615B5SOSM4601	441-0478
M	0615B5SOSM2602	441-0353	0615B5SOSM0902	441-0403	0615B5SOSM4602	441-0452
M/L	0615B5SOSM2603	441-0361	0615B5SOSM0903	441-0411	0615B5SOSM4603	441-0460
L	0615B5SOSM2604	441-0346	0615B5SOSM0904	441-0395	0615B5SOSM4604	441-0445
XL	0615B5SOSM2605	441-0387	0615B5SOSM0905	441-0437	0615B5SOSM4605	441-0486

PRO RANGE – OPEN TOE WITH FOOT

BE YOU Tonic Curvy Pro - 15-21mmhg						
	Blu Tonic		Nero		Ribes	
Size	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code
1 S-XL	0616B5SOSM2601	441-0734	0616B5SOSM0901	441-0775	0616B5SOSM4601	441-0817
2 M-XL	0616B5SOSM2602	441-0742	0616B5SOSM0902	441-0783	0616B5SOSM4602	441-0825
3 ML-XL	0616B5SOSM2603	441-0759	0616B5SOSM0903	441-0791	0616B5SOSM4603	441-0833
4 L-XL	0616B5SOSM2604	441-0767	0616B5SOSM0904	441-0809	0616B5SOSM4604	441-0841

CAPRI RANGE (FINISH BELOW THE KNEE)

BE YOU Capri - 15-21mmhg								
	Blu Tonic		Champagne		Nero		Ribes	
Size	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code
S	0621B5SOSM2601	440-8746	0621B5SOSM2201	440-8795	0621B5SOSM0901	440-8845	0621B5SOSM4601	440-8894
M	0621B5SOSM2602	440-8720	0621B5SOSM2202	440-8779	0621B5SOSM0902	440-8829	0621B5SOSM4602	440-8878
M/L	0621B5SOSM2603	440-8738	0621B5SOSM2203	440-8787	0621B5SOSM0903	440-8837	0621B5SOSM4603	440-8886
L	0621B5SOSM2604	440-8712	0621B5SOSM2204	440-8761	0621B5SOSM0904	440-8811	0621B5SOSM4604	440-8860
XL	0621B5SOSM2605	440-8753	0621B5SOSM2205	440-8803	0621B5SOSM0905	440-8852	0621B5SOSM4605	440-8902

BE YOU Capri Curvy - 15-21mmhg								
	Blu Tonic		Champagne		Nero		Ribes	
Size	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code	Product Code	PIP Code
1 S-XL	0622B5SOSM2601	440-9876	0622B5SOSM2201	440-9918	0622B5SOSM0901	440-9959	0622B5SOSM4601	440-9991
2 M-XL	0622B5SOSM2602	440-9884	0622B5SOSM2202	440-9926	0622B5SOSM0902	440-9967	0622B5SOSM4602	441-0007
3 ML-XL	0622B5SOSM2603	440-9892	0622B5SOSM2203	440-9934	0622B5SOSM0903	440-9975	0622B5SOSM4603	441-0015
4 L-XL	0622B5SOSM2604	440-9900	0622B5SOSM2204	440-9942	0622B5SOSM0904	440-9983	0622B5SOSM4604	441-0023

Use as part of a clinician-directed management plan for:

- Lower-limb lipoedema
- lipo-lymphoedema and Lower-limb oedema

Precautions and contraindications

Do not use in patients with acute infection/cellulitis, severe arterial insufficiency, uncontrolled cardiac failure, acute DVT, or known sensitivity to garment materials.

Use with caution and under clinical supervision in patients with fragile skin, open or highly exuding wounds, impaired healing, diabetes, or peripheral neuropathy.