

# credalast<sup>®</sup>

## Comfort

BRITISH STANDARD COMPRESSION HOSIERY

**Credalast<sup>®</sup> Comfort is an ideal compression therapy solution for the management of lower limb venous disease.**



**Versatile:** Credalast Comfort is available in Class 1 (14-17mmHg) and Class 2 (18-24mmHg) British Standard compression profiles. You are able to choose from below knee and thigh length with grip top in open and closed toe in 2 colours - black or beige. Our medical range of compression hosiery allows you to prescribe the appropriate level of support based on your patients' individual needs. Our range caters to a wide spectrum of medical conditions, ensuring comprehensive care for your patients.

**Ease of Use:** Credalast comfort was designed to make it easier for patients to put on and take off their hosiery. This means enhanced convenience for your patients, leading to improved compliance and overall satisfaction.


**Superior Fit:** The Lycra<sup>®</sup> Adaptiv yarn allows our hosiery to adapt to different limb shapes, providing a better fit experience. Patients will appreciate the support and enhanced comfort this unique yarn technology offers.

**Comfort in Motion:** Credalast Comfort ensures unrestricted movement and exceptional comfort throughout the day.

**Stay-in-Place Design:** We have increased the width of our top bands to reduce the need for constant adjustment. This includes the use of a NEW silicone beaded grip top.



### Measuring and selecting the correct size hosiery



Where possible, measurements should be taken early in the day or when the leg has been rested and swelling is at a minimum.

Please take each measurement from the specified points and match them against the sizing chart below:

**Below Knee = cC** (widest part of calf), **cB** (narrowest point above the ankle bone) and **L-Z** (foot length)

**Thigh Length = cG** (measure thigh circumference 5cm below gluteal fold), **cC** (widest part of calf), **cB** (narrowest point above the ankle bone) and **L-Z** (foot length)

Sizing	1	2	3	4	5
	Small	Medium	Large	Extra Large	Extra Extra Large
<b>cG</b> Thigh	40-52 cm	44-56 cm	48-60 cm	50-64 cm	
<b>cC</b> Calf	30-38 cm	32-40 cm	35-43 cm	37-46 cm	40-48 cm
<b>cB</b> Ankle	19-25 cm	21-27 cm	23-29 cm	24-32 cm	25-33 cm
<b>L-Z</b> Foot Length	20-24 cm	22-26 cm	24-28 cm	26-30 cm	26-32 cm

The coloured thread knitted in the stocking denotes sizing colour.



#### Washing Instructions

Machine wash using mild detergent. Do not iron. Do not tumble dry. Do not bleach. Do not dry clean or spot clean with chemical cleaning fluids. Dry flat.

**Daylong**  
DIRECT.

For further support please contact us on **Tel 0115 932 0144**  
**Fax 0115 944 0437** · Email [sales@daylongdirect.co.uk](mailto:sales@daylongdirect.co.uk)

# Class 1 Credalast Comfort



## Indications for use

Class I – Light Support (Compression at ankle 14-17mmHg)

Indications: Superficial or early varicose veins, varicose veins during pregnancy. Other indications; post operatively and for long distance travel to help prevent deep vein thrombosis (DVT) in normal risk groups.

Colour	Codes	Small	Medium	Large	X-Large	XX-Large
<b>Below Knee Closed Toe</b>						
Beige ●	Product Code	CR-CAD1BE2C	CR-CAD1BE3C	CR-CAD1BE4C	CR-CAD1BE5C	CR-CAD1BE6C
	PIP-Code	428-1184	428-1176	428-1150	428-1192	428-1200
	NHS Supply Chain	EGD87130	EGD87272	EGD87239	EGD87240	EGD87296
<b>Below Knee Open Toe</b>						
Beige ●	Product Code	CR-CAD1BE2O	CR-CAD1BE3O	CR-CAD1BE4O	CR-CAD1BE5O	CR-CAD1BE6O
	PIP-Code	428-1275	428-1267	428-1168	428-1283	428-1291
	NHS Supply Chain	EGD87295	EGD87113	EGD87214	EGD87131	EGD87192
<b>Thigh Length Grip Top Closed Toe</b>						
Beige ●	Product Code	CR-CAG1BE2C	CR-CAG1BE3C	CR-CAG1BE4C	CR-CAG1BE5C	CR-CAG1BE6C
	PIP-Code	428-1374	428-1366	428-1358	428-1382	428-1390
	NHS Supply Chain	EGD87171	EGD87274	EGD87292	EGD87216	EGD87084
<b>Thigh Length Grip Top Open Toe</b>						
Beige ●	Product Code	CR-CAG1BE2O	CR-CAG1BE3O	CR-CAG1BE4O	CR-CAG1BE5O	CR-CAG1BE6O
	PIP-Code	428-1473	428-1465	428-1457	428-1481	428-1499
	NHS Supply Chain	EGD87132	EGD87135	EGD87293	EGD87098	EGD87276
<b>Below Knee Closed Toe</b>						
Black ●	Product Code	CR-CAD1BL2C	CR-CAD1BL3C	CR-CAD1BL4C	CR-CAD1BL5C	CR-CAD1BL6C
	PIP-Code	428-1234	428-1226	428-1218	428-1242	428-1259
	NHS Supply Chain	EGD87065	EGD87068	EGD87150	EGD87212	EGD87093
<b>Below Knee Open Toe</b>						
Black ●	Product Code	CR-CAD1BL2O	CR-CAD1BL3O	CR-CAD1BL4O	CR-CAD1BL5O	CR-CAD1BL6O
	PIP-Code	428-1325	428-1317	428-1309	428-1333	428-1341
	NHS Supply Chain	EGD87095	EGD87067	EGD87081	EGD87066	EGD87269
<b>Thigh Length Grip Top Closed Toe</b>						
Black ●	Product Code	CR-CAG1BL2C	CR-CAG1BL3C	CR-CAG1BL4C	CR-CAG1BL5C	CR-CAG1BL6C
	PIP-Code	428-1424	428-1416	428-1408	428-1432	428-1440
	NHS Supply Chain	EGD87193	EGD87172	EGD87097	EGD87133	EGD87242
<b>Thigh Length Grip Top Open Toe</b>						
Black ●	Product Code	CR-CAG1BL2O	CR-CAG1BL3O	CR-CAG1BL4O	CR-CAG1BL5O	CR-CAG1BL6O
	PIP-Code	428-1523	428-1515	428-1507	428-1531	428-1549
	NHS Supply Chain	EGD87151	EGD87275	EGD87243	EGD87297	EGD87215

# Class 2 Credalast Comfort



## Indications for use

Class II – Medium Support (Compression at ankle 18-24mmHg)

Indications: Varicose veins of medium severity, venous ulcer treatment and prevention of occurrence and recurrence, mild oedema, varicose veins during pregnancy. Other indications; post operatively and for long distance travel to help prevent deep vein thrombosis (DVT) in higher risk groups.

Colour	Codes	Small	Medium	Large	X-Large	XX-Large
<b>Below Knee Closed Toe</b>						
Beige ●	Product Code	CR-CAD2BE2C	CR-CAD2BE3C	CR-CAD2BE4C	CR-CAD2BE5C	CR-CAD2BE6C
	PIP-Code	428-1572	428-1564	428-1556	428-1580	428-1598
	NHS Supply Chain	EGD87218	EGD87219	EGD87241	EGD87129	EGD87247
<b>Below Knee Open Toe</b>						
Beige ●	Product Code	CR-CAD2BE2O	CR-CAD2BE3O	CR-CAD2BE4O	CR-CAD2BE5O	CR-CAD2BE6O
	PIP-Code	428-1671	428-1663	428-1655	428-1689	428-1697
	NHS Supply Chain	EGD87194	EGD87174	EGD87072	EGD87071	EGD87070
<b>Thigh Length Grip Top Closed Toe</b>						
Beige ●	Product Code	CR-CAG2BE2C	CR-CAG2BE3C	CR-CAG2BE4C	CR-CAG2BE5C	CR-CAG2BE6C
	PIP-Code	428-1770	428-1762	428-1754	428-1788	428-1796
	NHS Supply Chain	EGD87056	EGD87248	EGD87073	EGD87143	EGD87282
<b>Thigh Length Grip Top Open Toe</b>						
Beige ●	Product Code	CR-CAG2BE2O	CR-CAG2BE3O	CR-CAG2BE4O	CR-CAG2BE5O	CR-CAG2BE6O
	PIP-Code	428-1879	428-1861	428-1853	428-1887	428-1895
	NHS Supply Chain	EGD87203	EGD87261	EGD87137	EGD87202	EGD87201
<b>Below Knee Closed Toe</b>						
Black ●	Product Code	CR-CAD2BL2C	CR-CAD2BL3C	CR-CAD2BL4C	CR-CAD2BL5C	CR-CAD2BL6C
	PIP-Code	428-1622	428-1614	428-1606	428-1630	428-1648
	NHS Supply Chain	EGD87116	EGD87197	EGD87157	EGD87074	EGD87249
<b>Below Knee Open Toe</b>						
Black ●	Product Code	CR-CAD2BL2O	CR-CAD2BL3O	CR-CAD2BL4O	CR-CAD2BL5O	CR-CAD2BL6O
	PIP-Code	428-1721	428-1713	428-1705	428-1739	428-1747
	NHS Supply Chain	EGD87115	EGD87047	EGD87223	EGD87220	EGD87051
<b>Thigh Length Grip Top Closed Toe</b>						
Black ●	Product Code	CR-CAG2BL2C	CR-CAG2BL3C	CR-CAG2BL4C	CR-CAG2BL5C	CR-CAG2BL6C
	PIP-Code	428-1820	428-1812	428-1804	428-1838	428-1846
	NHS Supply Chain	EGD87222	EGD87158	EGD87117	EGD87283	EGD87262
<b>Thigh Length Grip Top Open Toe</b>						
Black ●	Product Code	CR-CAG2BL2O	CR-CAG2BL3O	CR-CAG2BL4O	CR-CAG2BL5O	CR-CAG2BL6O
	PIP-Code	428-1929	428-1911	428-1903	428-1937	428-1945
	NHS Supply Chain	EGD87085	EGD87253	EGD87159	EGD87144	EGD87263